



## INTRODUCING A NEW STANDARD IN FEEDING JEJUNOSTOMY ACCESS

### THE CLINTEC® JEJUNO-CATH® KIT

- ▼ **NEEDLE/GUIDEWIRE TECHNIQUE**  
Subserosal, intramural catheter tunnel gives advantages of the Witzel Procedure without removing the bowel lumen.
- ▼ **EASIER TUBE REMOVAL**  
Patented internal fixation disc is made of soft silicone—gently collapses for easy removal without surgery. Internal fixation disc sutures to external disk to hold catheter securely—the tube cannot fall out inadvertently.
- ▼ **OPTIMAL 8 FR SIZE**  
Administer elemental and whole protein formulas conveniently.
- ▼ **HELPS PREVENT CATHETER BLOCKAGE**  
Removable Smart Port™ connector helps trap large particles and undissolved medications before they enter the catheter.

*The patented JeJuno-Cath silicone internal fixation disk folds into itself for easy removal.*



*JeJuno-Cath is inserted directly into the jejunum and secured by sutures through both internal and external fixation discs.*

FOR MORE INFORMATION, CONTACT YOUR LOCAL CLINTEC REPRESENTATIVE OR CALL

Sample from  
Bartlett's Peak Communications, Inc.  
630-766-9446

